


# TECHNICAL DESCRIPTION

---



## TeraSpike Companion

Advanced Low Noise Current Amplifier for THz Sensors

**THE** new TeraSpike Companion is the first commercially available low-noise amplifier specially designed and optimized for THz time-domain spectroscopy applications based on photoconductive detectors and emitters - including of course the full line of Proteemics TeraSpike (TS) microprobes. In contrast to conventional current amplifiers, the TS Companion offers important additional functions that are typically required for THz-TDS-related tasks, including e.g. Connections for Tx and Rx devices, an optical alignment mode, a switchable bias function for all inputs as well as a fully electronic control for system integration. A fine selection of amplification in steps of a fifth of a decade is included for optimal utilization of the maximum possible signal-to-noise ratio in any THz system. As required for operation with fast-scanning THz TDS modules (e.g. TeraFlash Pro, Toptica), it also offers a comfortable bandwidth reserve that ensures unadulterated THz signal detection up to the highest amplification range.


### Key Features

- Two input channels for Tx and Rx connection
- Operation mode switchable between optical alignment & THz detection
- Trans-impedance (gain) variable in  $0.2 \times 10^x$  steps from  $1 \times 10^6$  to  $8 \times 10^8$  V/A
- Bandwidth DC/400 Hz to 40 kHz
- Dynamic bias voltage
- Local and remote control options
- Wide operating supply range


### Applications

- Proteemics near-field probe operation including transceiver, emitter, or detector probes
- General Terahertz photoconductive antenna operation
- Photodiode and photomultiplier amplifier
- Multichannel current acquisition system
- Preamplifier for lock-ins, A/D converters, etc.

## Block Diagram




## Interface description




Input side view

Input1: SMA, isolated, jack female  
Input2: SMA, isolated, jack female  
Detector bias output: shield of the input1 SMA & input2 SMA



Output side view

Output: SMA, jack female



Power supply connector

Socket: Lemo series 0B, 3-pin fixed socket  
Plug: Matching parts  
FGG. 0B.303.CLAD52  
Pin 1: GND Pin 2: -12 V Pin 3: +12 V

## Specifications

| Model                                      | LNA-40k-8E8  |
|--|--|
| Trans-impedance gain [V/A] <sup>†</sup>    | $1 \times 10^6 \dots 8 \times 10^8$ , switchable in steps of $0.2 \times 10^x$ |
| -3-dB-Bandwidth [kHz] <sup>†</sup>         | 40, for all gain settings. Switchable to 0.4                                   |
| Amplifier output voltage [V]               | $\pm 12$   |
| Output impedance                           | Open   |
| Max. output current [mA]                   | $\pm 40$   |
| External bias supply for THz device inputs | $\pm 12$ V, max. 120 mA  |
| Indicator LED                              | On: Overload, Off: Normal operation  |
| Power supply voltage                       | $\pm 5$ V... $\pm 14$ V (max.), $\pm 12$ V (recommended)                       |
| Power supply current                       | Typ. $\pm 50$ mA   |

# TeraSpike Companion


Advanced Low Noise Current Amplifier for THz Sensors

|                             |           |
|-----------------------------|-----------|
| Weight [g]                  | 120       |
| Operating temperature range | 0... 70°C |

Test conditions:  $V_s = \pm 12\text{ V}$ ,  $T_A \approx \text{Room Temperature}$ , Load Impedance = 10 M $\Omega$

†Can be customized on request. For further information please contact us at [info@protemics.com](mailto:info@protemics.com)

## Digital Control



|                                      |   |
|--------------------------------------|---|
| Model                                | LNA-40k-8E8   |
| Control input voltage range LOW bit  | -12 V... 0.8 V, recommended: GND/0 V                      |
| Control input voltage range HIGH bit | 2.4 V... +12 V, recommended: 5 V                          |
| Control input current range          | $\pm 1\text{ }\mu\text{A}$ @ 0.8 V (LOW) ... 2.4 V (HIGH) |

| Pin    | Description                             |
|--------|---|
| Pin 1  | Signal output (connected to SMA)        |
| Pin 2  | Supply +12 V                            |
| Pin 3  | Digital control input: AC/DC            |
| Pin 4  | Digital control input: Gain             |
| Pin 5  | Digital control input: Gain             |
| Pin 6  | Digital control input: Controls Input1  |
| Pin 7  | NC                                      |
| Pin 8  | Digital control input: Controls Input2  |
| Pin 9  | GND                                     |
| Pin 10 | NC                                      |
| Pin 11 | Digital control input: Gain             |
| Pin 12 | Digital control input: Gain             |
| Pin 13 | Channel for supplying voltage to Input1 |
| Pin 14 | Digital control input: Input1/Input2    |
| Pin 15 | Channel for supplying voltage to Input2 |

## Remote Control Mode

In remote mode the device is controlled using a Sub-D 15-pin interface:

| Pin    | Function                                | Mode [Input/Output]<br>Address [LOW/HIGH]: Function |
|--------|---|---|
| Pin 1  | Signal output (connected to SMA)        | Output  |
| Pin 2  | Supply +12 V                            | Output  |
| Pin 3  | Digital control input: AC/DC            | Input: HIGH: AC ... LOW: DC                         |
| Pin 4  | Digital control input: Gain             | see gain selection                                  |
| Pin 5  | Digital control input: Gain             | see gain selection                                  |
| Pin 6  | Digital control input: Controls Input1  | Input: LOW: GND ... HIGH: Voltage                   |
| Pin 7  | NC                                      | -   |
| Pin 8  | Digital control input: Controls Input2  | Input: LOW: GND ... HIGH: Voltage                   |
| Pin 9  | GND                                     | Input   |
| Pin 10 | NC                                      | -   |
| Pin 11 | Digital control input: Gain             | see gain selection                                  |
| Pin 12 | Digital control input: Gain             | see gain selection                                  |
| Pin 13 | Channel for supplying voltage to Input1 | Input   |
| Pin 14 | Digital control input: Input1/Input2    | Input: HIGH: Input ... LOW: Input2                  |
| Pin 15 | Channel for supplying voltage to Input2 | Input   |

## Manual Control Mode

In manual mode the device is controlled using the manual control plug:

| Switch/<br>Terminal Block | Function                                       | Address/Designation |
|---------------------------|--|---------------------|
| Switch 1                  | Input1/Input2                                  | HIGH/LOW            |
| Switch 2                  | GND/Voltage, to input1                         | LOW/HIGH            |
| Switch 3                  | GND/Voltage, to input2                         | LOW/HIGH            |
| Switch 4                  | DC/AC  | LOW/HIGH            |
| Switch 5                  | NC   | -                   |
| Switch 6                  | NC   | -                   |
| Switch 7                  | A0 Gain  | see gain selection  |
| Switch 8                  | A1 Gain  | see gain selection  |
| Switch 9                  | A2 Gain  | see gain selection  |
| Switch 10                 | A3 Gain  | see gain selection  |
| Terminal Block            | Stable supply +12 V                            | +12 V               |
| Terminal Block            | Ground   | GND                 |
| Terminal Block            | Stable supply +1 V                             | +1 V                |
| Terminal Block            | Ground   | GND                 |
| Terminal Block            | Signal output (connected to SMA)               | Signal              |
| Terminal Block            | V2IN1, Channel for supplying voltage to Input1 | V2IN1               |
| Terminal Block            | V2IN2, Channel for supplying voltage to Input2 | V2IN2               |
| Terminal Block            | Ground   | GND                 |

# TeraSpike Companion

Advanced Low Noise Current Amplifier for THz Sensors


Gain Selection chart:

| Gain            | Pin/Address |          |           |          |
|-----------------|-------------|----------|-----------|----------|
| (V/A)           | Pin 11/A1   | Pin 4/A0 | Pin 12/A3 | Pin 5/A2 |
| $1 \times 10^6$ | LOW         | LOW      | LOW       | LOW      |
| $2 \times 10^6$ | LOW         | LOW      | LOW       | HIGH     |
| $4 \times 10^6$ | LOW         | LOW      | HIGH      | LOW      |
| $8 \times 10^6$ | LOW         | LOW      | HIGH      | HIGH     |
| $1 \times 10^7$ | LOW         | HIGH     | LOW       | LOW      |
| $2 \times 10^7$ | LOW         | HIGH     | LOW       | HIGH     |
| $4 \times 10^7$ | LOW         | HIGH     | HIGH      | LOW      |
| $8 \times 10^7$ | LOW         | HIGH     | HIGH      | HIGH     |
| $1 \times 10^8$ | HIGH        | LOW      | LOW       | LOW      |
| $2 \times 10^8$ | HIGH        | LOW      | LOW       | HIGH     |
| $4 \times 10^8$ | HIGH        | LOW      | HIGH      | LOW      |
| $8 \times 10^8$ | HIGH        | LOW      | HIGH      | HIGH     |

## Typical Application Schematics

Some typical configuration examples for using the TeraSpike Companion in combination with TeraSpike microprobes are depicted in the following. Pin / Switch settings that need to be modified are marked in bold in the tables below.

### Scheme 1: Probe Alignment



In this scheme, the probe is biased with an externally supplied voltage and the generated photo-current is amplified by the TeraSpike Companion. This scheme is typically used for the probe alignment of a probe-beam to the TeraSpike's photo-switch. The bias voltage is supplied to the probe by TeraSpike Companion's input connector outer shield (red channel for bias) via pin 13 (for input1) or 15 (for input2). For the TeraSpike Companion to operate in this scheme the following control address needs to be given at 15 pins Sub-D connector.

Control Address for TeraSpike Companion in remote control mode, where input1 is used for probe alignment

| Pin             | 1 | 2 | 3   | 4     | 5    | 6 | 7   | 8   | 9 | 10    | 11    | 12   | 13 | 14 | 15 |
|-----------------|---|---|-----|-------|------|---|-----|-----|---|-------|-------|------|----|----|----|
| Control Address | * | * | LOW | Gain* | HIGH | - | LOW | GND | - | Gain* | Vbias | HIGH |    |    | *  |


*\*for more details please refer to, the remote control mode table*

Switch status for TeraSpike Companion in manual control mode, where input1 is used for probe alignment

| Switch/Terminal | 1 | 2 | 3 | 4 | 5 | 6 | 7  | 8  | 9  | 10 | Terminal | V2IN1                   |
|-----------------|---|---|---|---|---|---|----|----|----|----|----------|-------------------------|
| Switch Status   | I | I | O | O | - |   | A0 | A1 | A2 | A3 |          | Supply +1 V externally* |

*\*for more details please refer to, the manual control mode table*

## Scheme 2: Terahertz measurement



In this scheme, the TeraSpike microprobes are operating as THz antenna and the resulting current signals are amplified by the TeraSpike Companion. No bias voltage to the probe is required. For the TeraSpike Companion to operate in this scheme the following control address needs to be given at 15 pins Sub-D connector.

**Control Address for TeraSpike Companion in remote control mode, where input1 is used for Terahertz signal**

| Pin  | 1 | 2 | 3    | 4     | 5 | 6   | 7 | 8   | 9   | 10 | 11    | 12 | 13 | 14   | 15 |
|--|---|---|------|-------|---|-----|---|-----|-----|----|-------|----|----|------|----|
| Control Address  | * | * | HIGH | Gain* |   | LOW | - | LOW | GND | -  | Gain* |    | *  | HIGH | *  |
| *for more details please refer to, the remote control mode table |   |   |      |       |   |     |   |     |     |    |       |    |    |      |    |


**Switch status for TeraSpike Companion in manual control mode, where input1 is used for Terahertz signal**

| Switch/Terminal | 1 | 2 | 3 | 4 | 5 | 6  | 7  | 8  | 9  | 10 | Terminal | V2IN1 |
|-----------------|---|---|---|---|---|----|----|----|----|----|----------|-------|
| Switch Status   | I | O | O | I | - | A0 | A1 | A2 | A3 |    |          | X     |

*\*for more details please refer to, the manual control mode table*

**NOTE:** If a user wants to use input2 of the TeraSpike Companion instead of input1, the necessary changes need to be done as given in the remote control mode table and manual control mode table depending on the mode of operation.

## Scheme 3: Transceiver operation



In this scheme, the input1 of TeraSpike Companion is biasing the probe for Terahertz generation, while the input2 is used for measuring the Terahertz signal. For the TeraSpike Companion to operate in this scheme the following control address needs to be given at 15 pins Sub-D connector.

# TeraSpike Companion

Advanced Low Noise Current Amplifier for THz Sensors

Control Address for TeraSpike Companion in remote control mode,  
where input1 is used for biasing the probe and input2 is for measuring the Terahertz signal

| Pin  | 1 | 2 | 3    | 4     | 5    | 6 | 7   | 8   | 9 | 10    | 11    | 12  | 13 | 14 | 15 |
|--|---|---|------|-------|------|---|-----|-----|---|-------|-------|-----|----|----|----|
| Control Address  | * | * | HIGH | Gain* | HIGH | - | LOW | GND | - | Gain* | Vbias | LOW | *  |    |    |
| *for more details please refer to, the remote control mode table |   |   |      |       |      |   |     |     |   |       |       |     |    |    |    |

Switch status for TeraSpike Companion in manual control mode,  
where input1 is used for biasing the probe and input2 is for measuring the Terahertz signal

| Switch/Terminal  | 1 | 2 | 3 | 4 | 5 | 6 | 7  | 8  | 9  | 10 | Terminal | V2IN1                   |
|--|---|---|---|---|---|---|----|----|----|----|----------|-------------------------|
| Switch Status  | O | I | O | I | - |   | A0 | A1 | A2 | A3 |          | Supply +1 V externally* |
| *for more details please refer to, the manual control mode table |   |   |   |   |   |   |    |    |    |    |          |                         |

Annotations: I: Logic HIGH, O: Logic LOW, X: Don't care


NOTE: For Transceivers probe alignment, please refer to Scheme1 of TeraSpike Companion.

Have more questions about TeraSpike Companion configurations? please feel free to contact us at [info@protemics.com](mailto:info@protemics.com)


## Typical Characteristics

At  $V_{S+} = +12\text{ V}$ ,  $V_{S-} = -12\text{ V}$ ,  $R_L = 10\text{ M}\Omega$ ,  $T_A \approx \text{Room Temperature}$ .

The below graph shows the AC characteristics of input 1:



The below graph shows the AC characteristics of input 2:







# TeraSpike Companion

Advanced Low Noise Current Amplifier for THz Sensors

## Dimensions



## Accessory



Manual Control Plug (MCP)

EU Konformitätserklärung  
EU Declaration of Conformity

1. Gerätetyp/Produkt  
*Apparatus model/Product* TeraSpike Companion
2. Name und Anschrift des Herstellers  
*Name and address of the manufacturer* Protemics GmbH  
Otto-Blumentahl Strasse 25  
52074 Aachen
3. Die alleinige Verantwortung für die Ausstellung dieser Konformitätserklärung trägt der Hersteller. *This declaration of conformity is issued under the sole responsibility of the manufacturer.*
4. Gegenstand der Erklärung  
*Object of the declaration* TeraSpike Companion
5. Der oben beschriebene Gegenstand der Erklärung erfüllt die einschlägigen Harmonisierungsrechtsvorschriften der Union.  
*The object of the declaration described above is in conformity with the relevant Union harmonisation legislation.*  
  
RICHTLINIE 2014/30/EU DES EUROPÄISCHEN PARLAMENTS UND DES RATES vom 26. Februar 2014 zur Harmonisierung der Rechtsvorschriften der Mitgliedstaaten über die elektromagnetische Verträglichkeit  
  
DIRECTIVE 2014/30/EU OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility
6. Angabe der einschlägigen harmonisierten Normen, die zugrunde gelegt wurden, einschließlich des Datums der Norm, oder Angabe anderer technischer Spezifikationen, für die die Konformität erklärt wird, einschließlich des Datums der Spezifikation:  
*References to the relevant harmonised standards used, including the date of the standard, or references to the other technical specifications, including the date of the specification, in relation to which conformity is declared:*  
  
DIN EN 55032:2016-02; VDE 0878-32:2016-02 Elektromagnetische Verträglichkeit von Multimediageräten und -einrichtungen - Anforderungen an die Störaussendung (CISPR 32:2015); Deutsche Fassung EN 55032:2015  
  
DIN EN 55035:2018-04; VDE 0878-35:2018-04 Elektromagnetische Verträglichkeit von Multimediageräten - Anforderungen zur Störfestigkeit (CISPR 35:2016, modifiziert); Deutsche Fassung EN 55035:2017
7. Nicht zutreffend.  
*No applicable.*
8. Zusatzangaben  
*Additional information* --

Unterzeichnet für und im Namen von:  
*Signed for and on behalf of:*

Protemics GmbH  
Otto-Blumentahl Strasse 25  
52074 Aachen

Ort und Datum der Ausstellung:  
*Place and date of issue*


Aachen, 11.7.2022

Name und Funktion  
*name, function*

  
Dr. Nagel, CEO

# TeraSpike Companion

Advanced Low Noise Current Amplifier for THz Sensors



## Contact Information

**Protemics GmbH**  
Otto-Blumenthal-Straße 25,  
52074 Aachen, NRW, Germany  
Phone: +49 0241 8867140  
Email: [info@protemics.com](mailto:info@protemics.com)  
<https://www.protemics.com/>

---

## Important Notice and Disclaimers

The information contained herein is believed to be reliable; however, Proteomics GmbH makes no warranties regarding the information contained herein and assumes no responsibility or liability whatsoever for the use of the information contained herein. All information contained herein is subject to change without notice. Customers should obtain and verify the latest relevant information before placing orders for Proteomics products. The information contained herein or any use of such information does not grant, explicitly or implicitly, to any party any patent rights, licenses, or any other intellectual property rights, whether with regard to such information itself or anything described by such information. **THIS INFORMATION DOES NOT CONSTITUTE A WARRANTY WITH RESPECT TO THE PRODUCTS DESCRIBED HEREIN, AND QORVO HEREBY DISCLAIMS ANY AND ALL WARRANTIES WITH RESPECT TO SUCH PRODUCTS WHETHER EXPRESS OR IMPLIED BY LAW, COURSE OF DEALING, COURSE OF PERFORMANCE, USAGE OF TRADE OR OTHERWISE, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.**